

WHAT ARE THE NEXT STEPS?

Cell Types

Collect data for a **wider range of cell types** applicable to cell therapies.

Design Changes

Incorporate **multiple channels** with **distinct geometries**. Enables more effective data collection.

Novel Syringe System

Aid the development of a patented syringe design aimed to greatly **increase cell therapy efficacy**.

TLDR *(Want a summary?)*

We are making injections safer for cells!

Low cell survival during injection causes therapies to fail. We have developed a microfluidic device to test exactly how much shear rate cells can handle.

This helps us find safer limits, leading to more effective and efficient injection systems. This will get more healthy cells where they need to go.

ACKNOWLEDGEMENTS

We would like to thank our supervisor Professor James Moore Jr., as well as Daniel Watson and Raphael Delattre for their guidance and support during our project.


WANT TO LEARN MORE?



You can find more information, CAD files and our final report by simply scanning the QR Code!

Contact Us

Department of Bioengineering
Imperial College London
South Kensington Campus
London SW7 2AZ, United Kingdom
Email: james.moore.jr@imperial.ac.uk



SMALL
CHANNELS

**BIG
IMPACT**

**Designing a microfluidic device
to effectively evaluate cell
viability under shear rates**

IMPERIAL

WHAT'S THE PROBLEM?



Cell therapies are a revolutionary way to deliver life saving treatment.

Then why does it have low success rates?

95% of cells end up dead after travelling through the needle.


95%

This is largely due to high shear rates - fast fluid layers sliding past each other, which can damage cells.

OVERALL AIM


Determine the critical shear rate of a given cell type

How did we approach this?



WHY MICROFLUIDIC DEVICES?

LESS CELL MATERIAL
FINER CONTROL
EFFECTIVE COST



Microfluidic devices allow precise flow control and are able to replicate injection conditions, while using a minimal number of cells.

Controlled Shear Delivered

Validated Research Platform

Cells exposed to a relevant shear rate showed 87.3% reduced viability compared to controls

WHAT HAVE WE ACHIEVED?


Our simulations accurately guided device design, exposing cells to precise shear rates.

Cell Viability Impact

Created a microfluidic device for future shear rate research on different cells types

FINAL DEVICE DESIGN

After many iterations based on experimental results, the design converged onto several essential features:



- Advanced flow control enables precise regulation of inlet flow rate, ensuring efficient use of resources.
- Narrow channel creates a region of elevated shear rate.
- Sudden channel widening decelerates and disperses cells for imaging.
- Allows for excess fluids and cell material to exit.